


# 3-in-1



3in1

3

가

3

,  
3

가

1

2

( , )

1 40 $\mu$ m

:

1


3in1

1

가

가

가



가 1

가 1

95%

2

2 5 $\mu$ m

:

5%

2

가

5 $\mu$ m

2  
5 $\mu$ m

3 0.01 $\mu$ m

:



3

99.99%

0.01 $\mu$ m

Donaldson Japan

DOP TEST

99.991%,

0.00017ppm w/w,

0.021kg f/cm<sup>2</sup>

가

0.45 $\mu$ m

가

3in1